

Anaerobic Culture Jar

Anaerobic Culture Jar 3.5 Lit.



Model: HOST/ANAR-2

Transparent Anaerobic Culture Jar of 3.5.lit. Capacity which is made of Polycarbonate. The lid of the jar is made of Aluminum pressure die Casting fitted with Pressure Valve, Vacuum Valve, Safety Valve & compound Pressure cum vacuum gauge, three finger aluminum clamp and sealing "O" ring.

This is recommended only where inside air of jar is removed thru vacuum pump and CO₂ is filled in the jar thru CO₂ gas cylinder to generate anaerobic condition in the jar.

Also supplied with plain aluminum pressure die casted lid which can be fitted when readymade gas packs are used to generate the anaerobic conditions inside the jar.

When readymade gas packs are used, there is no need of lid fitted with dial gauge etc., for extra life of equipment.