

ANAEROBIC CULTURE JAR (SS):



ANAEROBIC CULTURE JAR B.T.L. TYPE Modified mcintosh and field patterns, all metal construction consisting of a stainless steel jar with top flange recessed for a toroidal rubber sealing ring and a cover held in place by a strong bridge clamp.

The cover is fitted with two needle valves for evacuation and a sachet containing a catalyst.

A side arm on the jar is provided for connecting to semi solid Lucas indicator. (The jar is supplied without semi solid Lucas indicator).

With Pressure cum Vacuum Gauge.